



Performance & Lifestyle



## Force Line

**New Premium  
Elastic Knitted Supports**

- Smart Design
- Optimal Compression  
& Support
- Anatomical Fit
- Superior comfort

## Common Features & Benefits

- 3D contoured knit for effective compression
- Anatomically correct design and compression profile
- Wrinkle-free, non-slip fit for comfortable wear
- Breathable, hypoallergenic material



**German Engineering  
& Knitting Technology**

### MANUFORCE®

Compressive elastic knitted wrist support for active stabilisation of wrist joint

### EPIFORCE®

Compressive elastic knitted elbow support with epicondyle silicone inserts

### GENUFORCE®

Knitted knee sleeve with intelligent compression profile for optimal support and comfort

### MALLEOFORCE®

Knitted ankle sleeve with anatomical profile and silicone inserts for malleoli protection

### ACHILLOFORCE AIR™


Elastic knitted ankle support with pneumatic aircells providing intermittent compression on the Achilles Tendon



## Malleoforce®

### INDICATIONS

- Post-operative or post-traumatic use
- Chronic soft tissue inflammation
- Chronic instability
- Sprains, Moderate ligament injuries
- Recurrent swelling, load-related pain in the foot and ankle
- Mild osteoarthritis



### Smart design

Correct ankle flexion angle is integrated in the knit design for a perfect fit

### Protection & comfort

Anatomically shaped silicone inserts protect the malleoli & provide a massaging effect



### Focused flexibility

The softness of the yarn and the knitting pattern do not hinder ankle flexion

### Anatomical shape

The heel profile provides support and correct fit

## Manuforce®

### INDICATIONS

- Chronic, post-traumatic or postoperative inflammation of wrist joint and/or soft tissues
- Wrist sprains, Tendonitis (tenosynovitis)
- Carpal Tunnel Syndrome
- Mild osteoarthritis, Rheumatic wrist pain



### Comfortable silicone insert

Reduces pressure on bony prominence and nerves



### Malleable palmar aluminum stay for stabilization

Can be removed following treatment of acute phase



### Removable strap

To adjust compression level

## Genuforce®

### INDICATIONS

- Post-operative or post-traumatic use
- Chronic soft tissue inflammation
- Recurrent swelling or load-related pain
- Chronic instability, mild ligament injuries, meniscus injuries



### Brace stability

Silicon lined stays ensure non-slip grip

### Popliteal comfort

The softness of the yarn and the knitted pattern of the popliteal area ensure comfort and wrinkle-free movement




### Patellar stabilisation

Flexible silicone buttress maintains the patella in position and provides a massaging effect

## Achilloforce Air™

### INDICATIONS

- Chronic or post-traumatic inflammation of the Achilles tendon area
- Achilles tendonitis
- Plantar Fasciitis



### Interconnected pneumatic aircells for intermittent compression

The aircells provide intermittent compression and a massaging effect on the Achilles tendon area while walking

### Removable elastic strap

To adjust compression level



### Optimal effect

Best worn in a sports shoe

## Epiforce®

### INDICATIONS

- Chronic, post-traumatic or postoperative inflammation of the elbow joint
- Tennis elbow (lateral epicondylitis)
- Golfer's elbow (medial epicondylitis)
- Mild osteoarthritis, Rheumatoid arthritis
- Bursitis



### Removable strap

To adjust compression level

### Comfort zone in elbow crease


For comfort of wear

### Silicone inserts


Help relieve pressure on the elbow muscles

## ORDERING INFORMATION


DONJOY® GENUFORCE®					
Part Number	Size	Thigh Circumference 15 cm above patella (A)		Calf Circumference 5 cm under patella (B)	
82-0027-1	1 = XS	25 – 28 cm	9¾ – 11"	37 – 40 cm	14½ – 15¾"
82-0027-2	2 = S	28 – 30 cm	11 – 11¾"	40 – 43 cm	15¾ – 17"
82-0027-3	3 = M	31 – 34 cm	11¾ – 13½"	43 – 46 cm	17 – 18"
82-0027-4	4 = L	34 – 37 cm	13½ – 14½"	46 – 49 cm	18 – 19½"
82-0027-5	5 = XL	37 – 40 cm	14½ – 15¾"	49 – 52 cm	19¼ – 20½"
82-0027-6	6 = XXL	40 – 43 cm	15¾ – 17"	52 – 55 cm	20½ – 21¾"
82-0027-7	7 = XXXL	43 – 46 cm	17 – 18"	55 – 58 cm	21¾ – 22¾"




DONJOY® MALLEOFORCE®			
Part Number	Size	Ankle Circumference	
82-0028-1	1 = XS	17 – 19 cm	6¾ – 7½"
82-0028-2	2 = S	19 – 21 cm	7½ – 8¼"
82-0028-3	3 = M	21 – 23 cm	8¼ – 9"
82-0028-4	4 = L	23 – 25 cm	9 – 9¾"
82-0028-5	5 = XL	25 – 27 cm	9¾ – 10½"
82-0028-6	6 = XXL	27 – 29 cm	10½ – 11½"



DONJOY® ACHILLOFORCE AIR™			
Part Number	Size	Ankle Circumference	
82-0252-1	1 = XS	17 – 19 cm	6¾ – 7½"
82-0252-2	2 = S	19 – 21 cm	7½ – 8¼"
82-0252-3	3 = M	21 – 23 cm	8¼ – 9"
82-0252-4	4 = L	23 – 25 cm	9 – 9¾"
82-0252-5	5 = XL	25 – 27 cm	9¾ – 10½"
82-0252-6	6 = XXL	27 – 29 cm	10½ – 11½"



DONJOY® EPIFORCE®			
Part Number	Size	Forearm Circumference	
82-0253-1	1 = XS	19 – 22 cm	7½ – 8¾"
82-0253-2	2 = S	22 – 25 cm	8¾ – 9¾"
82-0253-3	3 = M	25 – 28 cm	9¾ – 11"
82-0253-4	4 = L	28 – 31 cm	11 – 12¼"
82-0253-5	5 = XL	31 – 34 cm	12¼ – 13½"
82-0253-6	6 = XXL	34 – 37 cm	13½ – 14½"



DONJOY® MANUFORCE®				
Part Number RIGHT	Part Number LEFT	Size	Wrist Circumference	
82-0254-1	82-0255-1	1 = XS	13 – 14.5 cm	5 – 5¾"
82-0254-2	82-0255-2	2 = S	14.5 – 16 cm	5¾ – 6¼"
82-0254-3	82-0255-3	3 = M	16 – 17.5 cm	6¼ – 7"
82-0254-4	82-0255-4	4 = L	17.5 – 19 cm	7 – 7½"
82-0254-5	82-0255-5	5 = XL	19 – 20.5 cm	7½ – 8"
82-0254-6	82-0255-6	6 = XXL	20.5 – 22 cm	8 – 8¾"



DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM

[www.DJOglobal.eu](http://www.DJOglobal.eu)



<sup>+</sup>  
MOTION IS MEDICINE™